# Donegal HA 2.0

40 mg/2 ml

Acido ialuronico sale sodico 2,0%  
2,0% hyaluronic acid sodium salt  
Hyaluronsäure Natriumsalz 2,0%

Revisiione: 12/2022

## Indicazioni

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## Indications

DONEGAL HA 2.0 acts as a substitute of synovial fluid into joints affected by degenerative or mechanical osteoarthritis that gives pain or reduced mobility.

## Description

DONEGAL HA 2.0 is a CE-marked medical device conforming to directive MDD 93/42/CEE, containing 20 mg/ml of hyaluronic acid obtained by fermentation and chemically unmodified. DONEGAL HA 2.0 is a transparent, sterile, non-pyrogenic and viscoelastic solution supplied in a 2 ml syringe. Hyaluronic acid is a natural polysaccharide present in numerous human tissues, especially synovial fluid, and acts on joints both as a cartilage and ligament lubricant and as a shock absorber. Many studies have shown that the injection of hyaluronic acid into joints affected by osteoarthritis restores synovial fluid viscosity and elasticity, with a consequent reduction in pain and improvement in the joint's mobility. DONEGAL HA 2.0 acts only in the joint into which it is injected, without exercising any systemic action.

## Composition

Main ingredient: hyaluronic acid sodium salt 2.0%. Other ingredients: sodium chloride, sodium phosphate and water for injectable preparations (to volume).

## Pack

DONEGAL HA 2.0 packages with 1 and 3 glass syringes each containing 40 mg of hyaluronic acid in 2 ml of saline solution buffered with sodium chloride in one blister foil. The content of the syringe is sterilized with moist heat.

## Instructions for use

Remove any fluid build-up in the joint before injecting DONEGAL HA 2.0, using the same needle for both removal and injection. Remove the protective cap from the syringe, taking care to avoid contact with the opening. Firmly screw a needle of 18-22 G diameter into the Luer-type lock neck to ensure a watertight seal. Before injection, treat the site with a suitable antiseptic solution. Inject DONEGAL HA 2.0 using an aseptic technique. Inject into the joint cavity only.

## Administration

DONEGAL HA 2.0 must be administered weekly for a total of 3 weeks or in any case according to the medical prescription.

## Warnings

- The pre-filled syringe contents are sterile. The syringe is packaged in a sealed blister pack. The outer surface of the syringe is not sterile.
- Do not use DONEGAL HA 2.0 after the expiry date reported on the pack.
- Do not use DONEGAL HA 2.0 if the pack or syringe is opened or damaged.
- The injection site must be in healthy skin.
- Do not inject into the bloodstream. Do not inject outside of the joint cavity, into tissue or in the synovial sacs.
- DONEGAL HA 2.0 has not been tested in pregnant or breastfeeding women.

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#### Distributore

Chiesi Italia S.p.A. Via Giacomo Chiesi, 1 43122 - Parma (Italia)

- DONEGAL HA 2.0 is single use and must not be re-sterilised.
- Avoid the simultaneous administration of DONEGAL HA 2.0 with other products for intra-articular use, in order to prevent any interaction.
- Do not administer DONEGAL HA 2.0 in the presence of copious intra-articular fluid build-up.
- Once the pack is opened, DONEGAL HA 2.0 must be used immediately and disposed of after use in conformance with applicable laws.
- Keep away from children.
- As with all invasive joint treatments, the patient is advised to avoid all intense physical activity for the first 2-3 days after the injection.

#### Contraindications

DONEGAL HA 2.0 must not be administered:

- to patients with known sensitivity to hyaluronic acid and its compounds
- in the case of skin infection or disease near the injection site
- if the joint is infected or highly inflamed.

#### Side effects

Various transitory reactions may follow DONEGAL HA 2.0 injection, such as pain, rigidity, sensation of heat, reddening or swelling. These secondary manifestations may be relieved by the application of ice to the treated joint. They normally disappear after a short time. If these symptoms persist, contact a doctor.

As for any invasive joint treatment, septic arthritis may arise if the necessary precautions are not observed during injection or if the injection site is not aseptic.

#### Storage and expiry

DONEGAL HA 2.0 must be stored in its original packaging at a temperature between 0 and 25°C, away from direct sunlight and frost. The expiry date is reported on the pack.

INTRA-ARTICULAR INJECTION MUST BE PERFORMED SOLELY BY A DOCTOR OR IN CONFORMANCE WITH LOCAL LEGISLATION.

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